


# BestScan®

## S3 IP67 Wireless Probe + S3 IP68 Pad

- Submersability
- Durability
- Longevity



※IP68 certified - dust tight and water proof (up to 1m for 30 mins)



### Applications

## BestScan® S3 Wireless

#### Main unit:

Sector probe, multifrequency 2.5 / 3.5 / 5.0 Mhz

Rechargeable Battery inside

AC Adapter for systems AC110V-250V 50/60H

Aluminum allow carry case

24 months for the main unit 12 months for probes

#### Optional:

Pad-7.0 inch Android Pad tablet. High resolution, special quick screenshot saving images

Pad-7.0 inch Android tablet. Water-proof dustproof corrosion-resisting IP68

**I**<sup>2</sup>  
ntelligent  
nsemination

**KUBUS**  
We Are Animal Reproduction

## Features of BestScan® S3 Wireless Probe

### Best-Scan® S3 Wireless broadband curved transducer

- 5 to 2 MHz extended operating frequency range
  - High resolution imaging for deeper applications
- Best-Scan® S3 Wireless works with your compatible smart device, giving you an incredible ultrasound system on top of the device's full functionality.

### Incredible image quality

Our portable ultrasound transducers and app include decades of expertise and innovation in ultrasound imaging to help you make fast, informed decisions.

*BMV's ultrasound app brings diagnostic capabilities to your compatible smartphone or handheld device.*



Model: S3 IP67 Wireless Probe

2.5 hours working time

Model: S3 IP68 Pad

 7000mAh lithium battery

 7.0 inch display



Shockproof



Anti-Corrosion



Compressive Force



Water Resistant



Dust-Proof



Anti-Scratch

The versatile ultrasound with wireless probe is made for pig-sheep-goat animal scanning

KUBUS

C/Varsovia, 20. Las Rozas de Madrid. 28232. España

916360268 [www.kubus-sa.com](http://www.kubus-sa.com) [kubus@kubus-sa.com](mailto:kubus@kubus-sa.com)



**I**<sup>2</sup>ntelligent  
nsemination

**KUBUS**  
We Are Animal Reproduction